



Interchim® Auto-Sampler puriFlash®AS-1 *Preparative LC Auto-Sampler*

Through puriFlash® system:
Injection
Liquid injection
w/column equilibration.




Load & Go™ technology:
through automated
6 port electric valve

Injection volume:
from 500µL up to 300mL

Needle for septum piercing
- 3 external solvent ways
- sample injection and wash
- bracketing injection mode
- Needle washing

Temperature control of the rack:
4°C to RT (option)

1 rack /slots for different test tubes capacity:
96 deep well plates, 2.5mL vials....
up to 1 liter bottles & custom racks



W: 14" - 37.0cm
D: 15" - 39.0cm
H: 18" 46.0cm

P/N: LO8850

